

# EUROCOR

---

CARDIOVASCULAR AND ENDOVASCULAR  
TREATMENT DEVICES

## **DIOR® DRUG-ELUTING PTCA BALLOON TECHNOLOGY SPECIFICALLY DESIGNED FOR CORONARY INTERVENTIONS**

Paclitaxel Paclitaxel is an active ingredient that inhibits cell replication; thus blocking the breakdown of microtubules during the metaphase and anaphase phases of mitosis. By selectively inhibiting the proliferation of smooth muscle cells, paclitaxel does not affect non-proliferating cells.

## **FREEWAY™ DEB – DRUG ELUTING BALLOON FOR PERIPHERAL SOLUTIONS – ONE STEP AHEAD**

Broad DEB OTW product portfolio available for de-novo, restenotic and AV fistula lesions:

**FREEWAY™ 014** – Designed specifically for the treatment of subgenual surgeries below the knee 0.014" OTW – Excellent balloon and crossover profile, smooth tip, very fast deflation time.

**FREEWAY™ 035** – Designed specifically for the treatment of femoropopliteal arteries 0.035" OTW – Excellent balloon with good pushability and secondary profile – Amorphous, non-crystalline bioshell coating Balloons 2–8 mm – Balloon lengths 20 –230 mm antegrade and crossed version.

